vitalmaxx®

BodiPro Massagegerät

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DE

BodiPro Massage Device

Instruction manual (Page 18)

GB

Appareil de massage BodiPro

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FR

BodiPro Massageapparaat

Handleiding (Pagina 48)

NL

Z 01467_VI

Inhalt

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Contents

Meaning of the Symbols in These Instructions _____	19
Safety Instructions _____	20
Intended Use _____	20
Warranty Terms _____	24
Items Supplied _____	24
Device Overview _____	25
Putting the Device into Operation _____	25
Unpacking _____	25
Inserting Batteries _____	26
Attaching / Removing the Attachments _____	27
Switching the Device On and Off _____	27
Use _____	27
Body Massage _____	28
Facial Massage _____	29
Cleaning and Storage _____	30
Troubleshooting _____	31
Technical Data _____	31
Disposal _____	32

Dear Customer,

We are delighted that you have chosen the “vitalmaxx” BodiPro Massage Device. With this battery-powered device, you can comfortably massage your body and face.

We hope you have a lot of fun with your new “vitalmaxx” BodiPro Massage Device.

If you have any questions, contact the customer service department via our website: www.service-shopping.de

Before using the device for the first time, please read the instructions through carefully and store them in a safe place. If the device is given to someone else, it should always be accompanied by these instructions. The manufacturer and importer will not accept liability if the directions in these instructions are not observed!

We reserve the right to modify the product, packaging or enclosed documentation at any time in conjunction with our policy of continuous development.

Meaning of the Symbols in These Instructions

All safety instructions are marked with this symbol. Read these instructions carefully and follow them to avoid injury to persons or damage to property.

Tips and recommendations are marked with this symbol.

Safety Instructions

Intended Use

- ☐ This device is intended to be used to massage the face and body. It is NOT suitable for therapeutic purposes.
- ☐ Only ever use the device directly on the skin and not on clothing. This could become caught in the device.
- ☐ This device should not be used by people (including children) with impaired sensory or mental faculties or who do not have adequate experience and/or knowledge of such devices, unless they are supervised by a person responsible for their safety, or have received instructions from such a person on how to use the device. Children must be supervised in order to ensure that they do not play with the device.
- ☐ Do not use the device while driving a car.
- ☐ This device is for personal use only and is not intended for commercial applications.
- ☐ Use the device only as described in the operating instructions. Any other use is deemed to be improper.

Danger of Injury

- ☐ Do not use the device in locations with high humidity, e.g. a bathroom. Do not use the device while you are showering or in the bath. There is the risk of an electric shock.
- ☐ Do not use the device if you have some form of electronic transplant device such as a pacemaker.
- ☐ If you have any health concerns or belong to a risk group, consult your doctor before using the device.

- ☐ Do not use the device if you are connected to some form of electronic medical equipment such as an ECG machine.
- ☐ The device must not be used by people who are connected up to some form of medical life support machine.
- ☐ People with heart complaints, circulation problems, varicose veins, diabetes, epilepsy, cancer, nerve damage or paralysis may only use the device following consultation with a doctor.
- ☐ Pregnant women must not use the device on their stomach.
- ☐ The device must not be used on children or on people who are asleep or unconscious. You must remain awake when using the device.
- ☐ People who have no feeling on their skin must not use the device.
- ☐ Do not use the device if you have a fever.
- ☐ Do not use the device on injured, infected, loose or (sun-) burnt skin, or on skin which has undergone plastic surgery, or in the immediate vicinity of acne, moles or varicose veins. Do not use the device if you have pains or swellings, suffer from haemophilia, diabetes, an immunodeficiency disorder or rosacea.
- ☐ If you experience increased pain or increased tension, skin irritations or hypersensitivity, switch off the device and stop using it. Consult a doctor.
- ☐ Keep the device away from your head, mouth, eyes, neck, joints and genitals.
- ☐ Do not use the device directly over your heart or chest.

- ☐ Do not use the device on tender, thin skin! You could be injured.
- ☐ Before every use, check whether the device and the accessories are in good condition. The device must not be used if it is deformed or damaged.
- ☐ Do not press the device too firmly onto your skin. Only ever apply slight pressure.
- ☐ Make sure that the rotating parts of the device only come into contact with the skin to be treated. Do not touch it with your fingers!
- ☐ During use, keep the device away from both head hair, eyelashes, eyebrows and from clothing, brushes, wires, cables or shoelaces to avoid any risks of injury.
- ☐ For reasons of hygiene, an attachment should only be used by one person. We do not recommend that you share it with other people.
- ☐ Do not use the device together with other electrical equipment, e.g. electric blankets.

Danger of Fire

- ☐ Do not use the device if highly flammable gases are present in the air.

Danger from Batteries

- ☐ Batteries can be deadly if swallowed. Therefore keep batteries and the device out of the reach of small children. If a battery is swallowed, seek medical assistance immediately.

- ☐ Before inserting the batteries, check whether the contacts in the device and on the batteries are clean; clean them if necessary.
- ☐ Use only the battery type specified in the technical data.
- ☐ Always replace all the batteries at the same time. Do not use batteries of different types or brands, or batteries with different capacities. When replacing batteries, note the correct polarity (+/–).
- ☐ Take the batteries out of the device when they are discharged or if the device is not going to be used for a prolonged period. This will prevent damage that might occur if a battery leaks.
- ☐ Should a battery have leaked, do not allow your skin, eyes and mucous membranes to come into contact with the battery acid. If you touch battery acid, rinse the affected sites with abundant fresh water immediately and contact a doctor as soon as possible.
- ☐ Remove a dead battery from the device immediately. Clean the contacts before you insert a new battery.
- ☐ The batteries must not be recharged or reactivated by other means, dismantled, thrown into a fire, immersed in fluids, or short-circuited.

Beware of Material and Device Damage

- ☐ In the event of damage/faults, switch the device off immediately.
- ☐ In the event of any damage, please contact the customer service department or a specialist. Do not try to repair the device or to replace any individual parts yourself to ensure your own safety.

- ☐ Only use original accessory parts.
- ☐ Do not expose the device to extreme temperatures, strong temperature fluctuations, direct sunlight, dust or moisture.
- ☐ Keep the device a sufficient distance away from sources of heat such as hobs or ovens.
- ☐ Do not drop the device or allow it to be knocked violently.
- ☐ Ensure that the device cannot fall into water or other liquids and that the main unit cannot become wet.
- ☐ Keep the device in a dry place at temperatures of between 0 and 40 °C.
- ☐ Do not allow the device to operate unattended. Switch off the device when you are not using it.
- ☐ Do not turn an attachment fitted on the device by hand (with the device switched off).

Warranty Terms

The following are excluded from the warranty: all defects caused by improper handling, damage or unauthorised attempts to repair. The same applies for normal wear.

Items Supplied

- 1x main unit
- 1x attachment for facial massage
- 1x attachment for body massage
- 4x 1.5 V batteries, AA type

Device Overview

- 1 Main unit
- 2 Battery compartment
- 3 Attachment for facial massage
- 4 Attachment for body massage
- 5 On/Off switch

Putting the Device into Operation

Unpacking

- Remove the device and accessories from the packaging. Check the items supplied for completeness and transport damage. **If you find transport damage, do not use the device (!), but contact the customer service department immediately.**
- Remove any films and stickers that may be present.

Inserting Batteries

1. Remove the battery compartment cover from the main unit in the direction shown by the arrow.

2. Pull the battery holder up and out of the battery compartment and insert four 1.5 V batteries of the AA type into the holder. Pay attention to the correct polarity (+/-).

3. Insert the battery holder into the battery compartment in such a way that the two metal contacts of the holder rest against the contacts in the battery compartment. The holder must lock in place in the battery compartment.

4. Slide the battery compartment cover back onto the main unit. Press down slightly as you do this so that it locks in place.

Replace the batteries in good time when the power of the device starts to become weaker.

Attaching/Removing the Attachments

- Take the attachment you want and plug it as far as it will go onto the axis of the main unit. It must engage audibly.
- To remove the attachment again, simply pull it off the main unit.

Switching the Device On and Off

- Slide the On/Off switch to the **ON** position to switch on the device.
The mounted attachment will start to rotate.
- To switch off the device, slide the switch to the **OFF** position.

Use

CAUTION!

- ☐ **Do not use the device on your chest** as this could cause cardiac arrhythmias!
- ☐ Check the attachments before using them. They must be in perfect condition.
- ☐ Make sure that no hairs can become caught in the attachment. As body hairs can also become caught on the attachment, they must be removed prior to the massage.
- ☐ Rub cream into the skin prior to the massage.
- ☐ If any pain or skin irritations occur during use, stop using the device immediately!


- ❑ Do not use the device on one spot for more than 1 minute and do not press the device too firmly against the skin in order to prevent skin irritations or injuries.
 - ❑ Do not massage more than 6 body parts one after the other. You should then take a break.
-

Body Massage

1. Place the attachment for body massage onto the main unit.
2. Apply massage oil, cream or gel to the part of the skin which is to be massaged and then wipe your hands dry.
3. Switch on the device. Place the rotating attachment onto the area which is to be massaged and move the device over the skin with circular movements. Move the device slowly and with only slight pressure. Massage one spot for no more than 1 minute.

- Make large circular movements when massaging the waist and stomach. Leave out the belly button!

- Make small circular movements when massaging the feet, legs, hips and bottom. Move the device slowly upwards (from the feet to the hips or from the thigh to the bottom).
4. Switch off the device after the treatment.

5. Remove the attachment from the main unit and clean it with clear, warm water. Allow the attachment to dry completely before you put it away.
6. If necessary, remove any flakes of skin dislodged by the massage or excess massage oil from your skin with clear, warm water and mild soap.

Facial Massage

1. Place the attachment for facial massage onto the main unit.
2. Tie your hair together. Make sure that there is no hair hanging over your face which may get caught up in the device. Perform the facial massage in front of a mirror to make sure that you do not get too close to your hairline with the device.
3. Apply massage oil, cream or gel to the face and then wipe your hands dry.

4. Switch on the device. Place the rotating attachment onto the lower part of the face and move it carefully up over the cheeks. Massage each half of the face for no more than 1 minute.
5. Then switch off the device.
6. Remove the attachment from the main unit and clean it with clear, warm water. Allow the attachment to dry completely before you put it away.
7. If necessary, remove any excess massage oil from your skin with clear, warm water and mild soap.

Cleaning and Storage

CAUTION!

- ☐ Make sure the device is switched off.
- ☐ No liquid must be allowed to enter the main unit. Do not immerse it in water or clean it under running water!
- ☐ Do not use any caustic or abrasive cleaning agents or cleaning pads to clean the device. These could damage the surface.

For reasons of hygiene, the device and the attachments used should be cleaned after every use. To do this, remove the attachment which has been used from the main unit.

- If necessary, wipe down the main unit with a slightly damp, soft cloth. Then dry it thoroughly.
- Rinse the attachments with clear, warm water. Allow the attachments to dry thoroughly before you store them away in a dry, clean place that is out of the reach of children and animals.

Troubleshooting

Problem	Possible cause/solution
The device is not working.	<ul style="list-style-type: none">• Are the batteries too weak or incorrectly inserted? If necessary, insert new batteries (see section "Inserting Batteries").• Has something become caught in the device that is now blocking the attachment?
The device heats up.	<ul style="list-style-type: none">• Have you used the device for too long without stopping?• Have you applied too much pressure?• Has something become caught on the attachment? <p>Allow the device to cool down for 1 hour before you switch it on again.</p>

Technical Data

Article No.: Z 01467

Model No.: MPI00012

Voltage supply: 4x 1.5 V batteries, AA type

Disposal

The packaging can be recycled. Please dispose of it in an environmentally friendly manner and take it to a recycling centre.

Dispose of the device in an environmentally friendly manner. It should not be disposed of in the normal domestic waste. Dispose of it at a recycling centre for used electrical and electronic devices. You can obtain more information from your local authorities.

Non-rechargeable and rechargeable batteries should be removed prior to disposal of the device and disposed of separately from the device. To protect the environment, non-rechargeable and rechargeable batteries may not be disposed of with the normal household waste, but must be taken to suitable collection points. Also note the relevant statutory regulations governing the disposal of batteries.

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All rights reserved.
DS Produkte GmbH, Am Heisterbusch 1, 19258 Gallin, Duitsland